

Technology Backgrounder


Diagnostic X-ray

Radiography solutions with SkyFlow at this year's RSNA

- **SkyFlow** <http://www.usa.philips.com/healthcare/product/HCNOCTN339/skyflow>
- **MobileDiagnost wDR** - existing cleared release 2.0
<http://www.usa.philips.com/healthcare/product/HC712001/mobilediagnost-wdr>
MobileDiagnost wDR - new release with sliding column 510(k) pending
- **DigitalDiagnost** - release with motorized ceiling suspension
<http://www.usa.philips.com/healthcare/product/HC712220/digitaldiagnost-digital-radiography-system>

SkyFlow

Let your work flow with SkyFlow



With **SkyFlow**, Philips offers a novel technology that allows clinicians to combine the ease of a grid-less acquisition workflow with the contrast of a grid image for bedside chest radiography. It supports an efficient workflow by providing a consistent grid-like image impression for a wide range of patient types. Without any user interaction or implications on workflow, SkyFlow enhances image contrast as a grid would.

- Fast and smooth grid-less workflow with fully automatic image contrast enhancement
- Superb image quality with grid-like contrast for all patients types, including bariatric
- Outstanding patient focus with automatic operation, short exam times, and comfortable positioning

SkyFlow is available for Philips MobileDiagnost wDR Rel. 2 and DigitalDiagnost Rel. 4 systems.